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
connections Home Profiles Communities Apps Request a Community **Field Resources**

Service Assist

Wikis

- Welcome to Service Assist
- Popular Service Communities
- Popular Service Forums
- Popular Service Blogs

Index

Members

Trash

You are in: [Service Assist](#) > Welcome to Service Assist

Welcome to Service Assist

57 Like | Updated Jul 07 by Mary Wood | Tags: part_locator

Page Actions

Service Assist is the place for finding service knowledge.

- Service Assist
- File Sync
- Prod, MKT, & Field Ops Communities
- PIP Announcements
- Product Ordering Guides
- Site Prep Guides
- Product Brochures
- Specification Sheets

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■ Service Assist Wiki> Qualification/Compliance(检定/合规性)

Welcome to Service Assist

55 Like | Updated 7/7/17 by Mary Wood | Tags: part_locator

Page Actions ▾

Service Assist is *the* place for finding service knowledge.

Recently Released

Service Communications News

Popular

Key Documents

Part Locators & Reference Info

Installation Checklists

PM Protocols

PM Kits

Help

User Manuals

Care and Use Manuals

Site Prep Guides

Software Support

Troubleshooting

Video Library

Qualification/ Compliance

Contamination Control

Service Notes

Tools & Utilities

Spare Parts Search

Portal Search

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Qualification/Compliance			
The following is a list of Waters Qualification and Compliance information. This is a live document and will be updated as new information becomes available. If you have any questions or concerns regarding this information, contact Ginny Halfpenny at x8422 or Paul St. Jean at x2322.			
Qualifications Database	Compliance Services Lotus Notes Database	AQT Qualification	AQT Information on Waters.com
GxP Compliance	GxP Compliance on Waters.com		ACQUITY AQT Version 2.20: ACQUITY Maintenance Procedures ACQUITY OQ Tests Miscellaneous Forms
Parts/Kits	Compliance Kits - By Kit Description		Alliance AQT Version 3.00: 2475 MP 2487 MP 2695 MP 2795 MP 2996 and 996 MP RI Detector MP Miscellaneous Forms Release Notes, Protocols, and Test Examples Empower AQT SAT/IN AQT Version 3.00
SystemsQT Qualification	ACQUITY: ACQUITY Arc SystemsQT ACQUITY SystemsQT Documentation ACQUITY UPC² SystemsQT PATROL SystemsQT Alliance with Empower SystemsQT Empower AQT: Empower 1 Empower 2 Empower 3 Modular HPLC: Modular HPLC SystemsQT Documentation Empower 3 SQT Self-Directed SQT Using Empower SystemsQT Qualification Tool: ACQUITY UPLC Systems Agilent 1100/1200 Systems Alliance HPLC Systems	Qualification Workbooks	Qualification Workbook Information
		UNIFI SystemsQT	UNIFI SystemsQT Documentation
		Understanding Qualification	2487 TUV Wavelength Accuracy Qualification 2489 TUV Wavelength Accuracy Qualification 2996 Wavelength Accuracy Qualification

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 Qualification Workbooks			
The following is a list of Waters Qualification Workbooks and other qualification procedures and documents. This is a live document and will be updated as new information becomes available. If you have any questions or concerns regarding this information, contact Ginny Halfpenny at x8422.			
Informatics	Chromatography Data Software (CDS): Automated Development System Software (716001039 Rev D) Empower AQT (Listing of Empower AQT Service Notes) Empower 3 Enterprise/Workgroup for Servers (715002791 Rev B) Empower 2 Enterprise/Workgroup for Servers (71550066001 Rev A) Empower Enterprise/Workgroup for Servers (71550047001 Rev B) Empower Systems (71550032201 Rev B) Empower SQT for PATROL Millennium³² 4.0 Systems (71550011001 Rev A) Millennium³² 4.0 Systems for Servers (71550011812 Rev B) SAT/IN AQT 3.00 MassLynx: AppSight Black Box Software IO/OQ (71550064001 Rev A) MassLynx 4.0/4.1/4.2 Software Systems (71550093001 Rev E) MassLynx 4.1 w/Security NuGenesis: Empower+Archive Solution UNIFI: Regulated Bioanalysis System for UNIFI v1.5	MS Technology	ACQUITY QDa Detector: ACQUITY QDa w/ACQUITY Arc System (715004743 Rev B) ACQUITY QDa w/ACQUITY UPC² System (715004742 Rev B) Q-ToF: Q-ToF micro (71550002001 Rev B) Q-ToF Premier (71550012301 Rev B) Xevo G2 QToF (715002508 Rev D) Xevo G2-S QToF, MassLynx (715003605 Rev A) Xevo G2-XS QToF, MassLynx (715004593 Rev A) Xevo QToF (715002075 Rev B) Single Quadrupoles: 3100 Mass Detector SQ Detector 2, MassLynx (715003079 Rev A) SQ Detector 2, Empower (715005150 Rev A) ZQ 2000/4000, Empower/Millennium³² 4.0 (71550070012 Rev A) ZQ 2000/4000, MassLynx (71550024001 Rev B) SYNAPT: SYNAPT Family (71550230001 Rev A) SYNAPT G2 MS/HDMS, MassLynx (715002050 Rev C) SYNAPT G2-S, MassLynx (715003400 Rev A) SYNAPT G2-Si, MassLynx (715004241 Rev A) Time of Flight: GCT Premier, MassLynx (71550063001 Rev A) LCT Premier, MassLynx (71550008001 Rev A) LCT Premier XE, MassLynx (71550500001 Rev B) MALDI micro MX, MassLynx (71550009001 Rev A) Xevo G2 ToF, MassLynx (715002919 Rev B) Xevo G2-S ToF, MassLynx (715003605 Rev A) Xevo G2-XS ToF, MassLynx (715004593 Rev A) Tandem Quadrupoles: Quattro micro API, MassLynx (71550001001 Rev B) Quattro micro GC, MassLynx (71550019001 Rev A) Quattro Premier XE, MassLynx (71550505001 Rev A) Xevo TQ MS, MassLynx (715001928 Rev B) Xevo TQ-S, MassLynx (715002207 Rev D) Xevo TQ-S micro, MassLynx (715004603 Rev A)
LC Technology	ACQUITY UPLC: ACQUITY 30-cm Column Heater/Cooler (715002425 Rev A) ACQUITY AQT 2.20 ACQUITY UPLC w/FLR, MassLynx (71550220012/3/5 Rev A) ACQUITY UPLC w/TUV/PDA/2996, MassLynx (71550094001 Rev F) ACQUITY UPLC Family w/UV Detectors, MassLynx (715003394 Rev B) ACQUITY UPLC H-Class w/FLR, MassLynx (715002307 Rev B) ACQUITY UPLC H-Class w/TQD, UV, MassLynx (715002309 Rev A) ACQUITY UPLC H-Class w/UV, MassLynx (715002306 Rev A) ACQUITY UPLC I-Class w/UV, MassLynx (715003167 Rev B) LC Systems: Alliance 2690/2695 (71550269501 Rev E) Alliance 2695 (71550269512 Rev G) Alliance 2790/2795 Alliance 2796 (71550279601 Rev A) Alliance AQT 3.00 Alliance GPC 2000 (71550016901 Rev A)		

- IBM Connections (connect.waters.com) > Service Assist Wiki > Certificates of Structural Integrity(结构完整性证书)

Qualification/
Compliance

Contamination
Control

Service Notes

Tools & Utilities

Spare Parts
Search

Portal
Search

A to Z index Any other important service links to add? Contact Waters Customer Experience and Knowledge Management at service_assist@waters.com.

- Agilent
- Application Notes
- Care and Use Manuals
- Certificates of Structural Integrity
- Chemistry Returns
- CIT (Customer Interaction Tracking)
- Qualification/Compliance
- Qualification Workbooks (QWBs)
- Quick Install Guides
- Radiation Safety (EHS Manual)
- Reference Information
- Release Notes

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■ Waters.com > Support Library(支持文库)

The screenshot displays the Waters Support Library interface. At the top, the Waters logo and tagline 'THE SCIENCE OF WHAT'S POSSIBLE.®' are visible. Below this is a navigation bar with links for INDUSTRIES, PRODUCTS, RESOURCES, EDUCATION & EVENTS, SERVICES & SUPPORT, and ABOUT WATERS. A breadcrumb trail indicates the current location: Home > Services & Support > Support > Support Library. On the left, a sidebar menu lists various support categories, with 'Support Library' currently selected. The main content area is titled 'Support Library' and features a search bar with a 'Search' button. Below the search bar, it shows 'Find Documents: Content Type: Certificates of Structural Integrity' with a 'Clear All' option. A 'Product Area' dropdown menu is also present. The results section shows '1-25 of 443 Results' and a 'Go To' page selector set to 'Page 1'. Two document entries are listed under the 'Description' header: 'Headspace Control Option Version 3.0 Certificate of Structural Integrity' and 'ChromScope Versions 1.21, 1.41, and 1.5 Certificate of Structural Integrity'. Each entry includes its support number, publication date, content type, and a brief description of the certificate's coverage.

Waters
THE SCIENCE OF WHAT'S POSSIBLE.®

INDUSTRIES PRODUCTS RESOURCES EDUCATION & EVENTS SERVICES & SUPPORT ABOUT WATERS

Home > Services & Support > Support > Support Library

iRequest Technical Service /
Support
Support
Support Library
How To Videos
Spare Parts
Laboratory Analytics
Instrument Upgrades
Compliance Services
Instrument Services
Implementation Services
Software Services
Software Upgrades
Leasing and Financing Services
Service Excellence

Support Library

Search

Find Documents: Content Type: Certificates of Structural Integrity  **Clear All** 

Product Area 

1-25 of 443 Results Go To 

Description

[Headspace Control Option Version 3.0 Certificate of Structural Integrity](#)
Support Number: CSTR134840266 | Publication Date: April 10, 2015 | Content Type: Certificates of Structural Integrity
The attached Certificate of Structural Integrity covers Headspace Control Option Version 3.0 for Empower.

[ChromScope Versions 1.21, 1.41, and 1.5 Certificate of Structural Integrity](#)
Support Number: CSTR134837825 | Publication Date: March 31, 2015 | Content Type: Certificates of Structural Integrity
The attached Certificate of Structural Integrity covers ChromScope Software Versions 1.21, 1.41, and 1.5.

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- 05 Apr 2009
- 05 Apr 09
- 05 April 2009
- 05 April 09
- April 05, 2009

Certification		
The undersigned performer attests that the Waters ACQUITY UPLC System is installed according to the above procedure.		
Performer	<u>Stanoula Powell</u> Signature	Date <u>28 Oct 2013</u>

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Sample Organizer Heater/Cooler Temperature Accuracy Test Sample Organizer Not Installed SPP
Initials

Temperature Values Acceptable Range Test Results

Set-Point Temperature Observed Temperature +4.0 °C of Set-Point Temperature Pass Fail

Initials Initials Initials

Sample Organizer Plate Positioning Test

Position Number

Top Plate Position _____

Second Plate Position _____

Third Plate Position _____

Bottom Plate Position _____

N/A SPP 28 Oct 2013

To pass, each of the 4 plates must correctly transfer from the sample organizer, to the sample manager, and then back to their original position in the sample organizer.

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Waters Customer Assurance Documentation
US Field Tool Calibration Program
Informational

Rev: 02
Creation or revision date: 11/09/2012
Approval date: 11/14/2012

Originator: Jeff M Billingham

Model: Field tool calibration

The purpose of an Informational Service Note is to communicate detailed product related information to the Field Service Engineer.

Scope Technical Discussion

Below is a summary of the US Field **Tool Calibration** Program. It contains the contact information for our **Calibration** vendor and other pertinent information about the program. All FSE's are responsible for the timely recalibration of their **tools**. The attachment below, US **Tool** Cal Prgm_081712.docx, describes the US Field **Tool Calibration** Program. Also attached are two forms for use the the program.

Product Affected

Field **tool calibration**

For your information

Calibration Vendor Contact Information:

Tektronix Service Solutions
Attn: Natalie Tollberg
1946 Greenspring Drive
Suite A
Timonium, MD 21093

Account Administrator: Natalie Tollberg email: natalie.tollberg@tektronix.com phone: 410-842-1000
X1101

If you need expedited turn around, please contact Natalie. Return shipping is the standard FedEx 3 day unless you request otherwise.

Exchange Program is available for most **tools**.

Tektronix will send a Due for **Calibration** Recall Notice ~45 days prior to your next **calibration** due date and another reminder at 2 weeks prior to the due date if it has not been sent in for recalibration. The sender is "alerts" in your mailbox.

Our internal Waters worldwide field **tool calibration** process will independently generate Due for **Calibration** Notices that will be emailed to you at 90, 60 and 30 days prior to the **calibration** due date. If your **tools** become overdue for **calibration** they will disappear from your FSR dropdown list and you will get an Overdue for **Calibration** Notice at days 1, 20 and 40 past due. The sender is "Note Developer1" in your mailbox.

Service Action Required

All FSE's are responsible for the timely recalibration of their **tools**.

Implementation

N/A

Attachments

 [US Tool Cal Prgm_081712.docx](#)  [Exchange Pool Request Form_08Jan2013.docx](#)  [TSS Return to Depot Form.xls](#)

Here is a link to the Compliance Services page of the Waters portal. The Compliance Services page includes a link to the US field **tool calibration** process.
<http://portal.waters.com:57100/inj/cervlet/prtportal/prtroot/docca/guid/00bf82a1-7dac-2e10-b9ba-fc1ed93c7886>

THIS SECTION IS FOR THE ASIAN OPERATION AND IS TO BE USED FOR TRANSLATION INTO LOCATION LANGUAGE

THIS DOCUMENT AND ITS CONTENTS ARE EXCLUSIVE PROPERTY OF WATERS CORPORATION, MILFORD MASSACHUSETTS. THIS DOCUMENT MAY NOT BE REPRODUCED IN ANY FORM WHATSOEVER, WITHOUT PRIOR WRITTEN PERMISSION FROM WATERS CORPORATION.

This note is status is "Approved". The note was approved by Jeff M Billingham/Waters on 11/14/2012

Service Communications -- Use this section to request Service Communications to publish this information into a permanent document i.e: Service procedure, Tech note, installation checklist, part locator, or troubleshooting information

Please select the appropriate Contact person:

Service Communications Instrument Contact person ☐ Priscilla Lotfy/Waters

Service Communications Data - S/W Support Contact person ☐ Paul M Fiore/Waters

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[COMPLIANCE SERVICES FIELD GUIDELINES]

FSE SELF-PREPARATION STEPS

- ☐ Consult the Service Notes Database or your regional compliance services expert resource for the latest information on the services to be provided. (Note: If you don't have a regional expert, contact a Global Services Support Compliance expert at x8070, Option 6.)
- ☐ If there are applicable Validation Service Notes, print them for inclusion in the qualification record.
- ☐ Ensure that your training certificates and tool calibration certificates are current. Print them for inclusion in the qualification record.
- ☐ Confirm that the requisite materials have been ordered and review the timing and expected duration of planned engagement.

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

[COMPLIANCE SERVICES FIELD GUIDELINES]

FSE AND CUSTOMER DISCUSSION AND PREREQUISITES

Pre-Call or Onsite

- ☐ Discuss with the customer the purpose and scope of the testing.
- ☐ Discuss the details of the latest protocols and note customer preferences:
 - ☐ Systems to be qualified with the latest protocols (preferred).
 - ☐ Systems to be qualified with protocols currently in use at their site.
- ☐ Ensure the customer is alerted that pre-approval of the protocol is required. Offer to provide reference copies in advance.
- ☐ Confirm that you have your own account on the customer's system.
- ☐ Ensure the customer is informed that a project may be restored to their system during the Compliance Services engagement.
- ☐ Verify the customer is aware of required materials they are expected to provide including mobile phase, calibrated balances, pipettes, etc.

Onsite

- ☐ Confirm Good Documentation Practice requirements for date and time, signatures, and initials.
- ☐ Discuss where maintenance, repair, and qualification activities will be logged (for example, in the system maintenance procedures and log).

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

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Waters

THE SCIENCE OF WHAT'S POSSIBLE.™

Waters ConnectionsAQTSM for EmpowerTM 2 Software Qualification

License to Qualify

This document grants the purchaser the right to qualify **ONE** Empower 2 Software node (personal, enterprise, workgroup or LAC/E³²) in conjunction with the ConnectionsAQT for Empower 2 Software Advanced Qualification Technology software. There must be one license for each node that is qualified by the application.

Qualification of an Empower 2 Software node without a valid license constitutes a violation of this license.

To be filled in by Customer

This license is only valid for the Empower 2 Software node name listed below.

Node name: _____

Signature: _____ Date: _____
(Authorized Customer Signature)

Waters Corporation

Within the United States and Canada
Contact: Waters Technical Services: 1-800-252-4752

Outside the United States and Canada please contact your local Waters office

Please Retain this Document

P/N: 716001780 Rev A

Waters Empower SystemsQT for Software License to Qualify Empower 3 Software

This document grants the purchaser the right to qualify **ONE** (1) Empower 3 Software node (personal, enterprise, workgroup, LAC/E³² or server) using EmpowerTM SystemsQTSM Qualification Tool. There must be one license for each node that is qualified by the application. This license must be retained with the qualification documentation when printed reports are used or in a suitable repository when electronic reports are used.

Qualification of an Empower 3 Software node without a valid license constitutes a violation of this license. License duration should be indicated below at time of use and match that identified on the purchase order:

One Qualification ☐
Self Directed (one year) ☐
One Year Plan ☐
Two Year Plan ☐

This license is only valid for the Empower 3 Software node listed below.

Node name _____

Signature _____

Date _____

Within the United States and Canada

Contact: Waters Technical Services: 1-800-252-4752

Outside the United States and Canada please contact your local Waters office

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**Waters ConnectionsAQTSM for
Method Validation Manager Qualification
on EmpowerTM 2 Software**

License to Qualify

This document grants the purchaser the right to qualify **ONE** Method Validation Manager (MVM) Software node (on personal, enterprise, workgroup or LAC/E³²) in conjunction with the ConnectionsAQT for MVM Advanced Qualification Technology software. There must be one license for each node that is qualified by the application.

Qualification of a Method Validation Manager Software node without a valid license constitutes a violation of this license.

To be filled in by Customer

This license is only valid for the MVM Software node name listed below.

Node name: _____

Signature: _____ Date: _____
(Authorized Customer Signature)

Waters Corporation
Within the United States and Canada
Contact: Waters Technical Services: 1-800-252-4752
Outside the United States and Canada please contact your local Waters office

Please Retain this Document

P/N: 71600xxxx Rev A

**Waters EmpowerTM SystemsQTTM
SAT/IN Qualification**

License to Qualify

This document grants the purchaser the right to qualify **ONE** SAT/IN or eSAT/IN in conjunction with the Empower Systems Qualification Tool (SystemsQT) software. There must be one license for each SAT/IN or eSAT/IN that is qualified by the application. This license must be retained with the qualification documentation when printed reports are used or in a suitable repository when electronic reports are used.

Qualification of a SAT/IN without a valid license constitutes a violation of this license. License duration should be indicated below and match that identified on the purchase order:

One Qualification	☐
Self Directed (one year)	☐
One Year Plan	☐
Two Year Plan	☐

To be completed at time of use

This license is only valid for the SAT/IN or eSAT/IN Serial Number listed below.

SAT/IN or eSAT/IN Serial Number: _____

Date: _____

Please retain this document

Within the United States and Canada
Contact: Waters Technical Services: 1-800-252-4752
Outside the United States and Canada please contact your local Waters office

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

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Waters

THE SCIENCE OF WHAT'S POSSIBLE.™

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

LC/GC System Qualification

License to Qualify

This document grants the purchaser the right to **ONE** Liquid or Gas Chromatography System Qualification in conjunction with the Connections®AQT (Advanced Qualification Technology) software, or the Waters SystemsQT™ (Systems Qualification Tool). There must be one license for each system that is qualified by the application.

Qualification of a System without a valid license constitutes a violation of this license. This license is valid on date demand qualification service is performed or for the duration of an applicable service plan with qualification option.

Completed license to be signed by customer

This license is only valid for modules with the serial numbers listed below. Enter N/A if not applicable.

Solvent Manager or Separation Module Serial Number: _____

Sample Manager or Separation Module Serial Number: _____

Column Manager or Column Heater/Cooler Serial Number: _____

Sample Organizer Serial Number: _____

Detector Serial Number: _____

Signature _____ Date _____
(Authorized Customer Signature)

Waters Corporation
Within the United States and Canada
Contact: Waters Technical Services: 1-800-252-4752
Outside the United States and Canada please contact your local Waters office

Please retain this document in the
Qualification Documentation Binder

PN: 716002319 Rev A

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

Additional LC/GC Detector Qualification

License to Qualify

This document grants the purchaser the right to **ONE** additional Liquid or Gas Chromatography Detector Qualification in conjunction with the Connections®AQT (Advanced Qualification Technology) software for systems, or the Waters SystemsQT™ (Systems Qualification Tool). There must be one license for each system and additional detector that is qualified by the application.

Qualification of a system or detector without a valid license constitutes a violation of this license. This license is valid on date demand qualification service is performed or for the duration of an applicable service plan with qualification option.

Completed license to be signed by customer

This license is only valid for the detector module with serial number listed below.

Detector Serial Number: _____

Signature _____ Date _____
(Authorized Customer Signature)

Waters Corporation
Within the United States and Canada
Contact: Waters Technical Services: 1-800-252-4752
Outside the United States and Canada please contact your local Waters office

Please retain this document in the
Qualification Documentation Binder

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Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

**LC System Qualification
Service Certificate**

This document covers the labor for ONE Waters® Liquid Chromatography System Qualification in conjunction with the Connections®AQT (Advanced Qualification Technology) software and/or the Waters SystemsQT™ (Systems Qualification Tool).

Includes: Qualification Documentation, Data System Qualification, Organizational Binder and Qualification Overview documentation, Qualification Test Chemicals and materials as required, Qualification labor provided by a Waters Certified Compliance Specialist and Local (Zone 1) travel charges.

Pre-use Qualification Services for New Systems
Installation Qualification Documentation will be completed
Operational Qualification testing and documentation will be completed
System Performance testing will be completed and documented for systems with LC Detector

Post-use Qualification Services for Installed or Existing Systems
Maintenance documentation completed for Performance Maintenance service if purchased
Operational Qualification testing and documentation will be completed
System Performance testing will be completed and documented for systems with LC detector

This certificate is for documentation purposes only and cannot be used as credit towards any other product or service.

741000348

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

**LC System Module
Qualification Labor Certificate**

This document covers the Qualification Service labor for ONE Waters® Liquid Chromatography System Module

Includes: Qualification labor provided by a Waters Certified Compliance Specialist, and Local (Zone 1) travel charges

This Certificate does not cover:
Time and materials necessary to prepare an instrument for qualification.
Workbooks or protocols used during qualification.
Calibration test kits maintenance kits, repair parts, chemicals, mobile phase or column required for qualification.
Additional travel expenses for a customer site that is outside of Waters Zone 1 travel regions.

Additional Information
All qualification workbooks, protocols, calibration test solutions and columns should be purchased by the customer.
These materials should be on-site before the qualification service is performed.
Additional travel expenses and service materials (if required) will be billed at current rates during the qualification visit.
Mobile phase and proper laboratory facilities are the responsibility of the customer.

This certificate is for documentation purposes only and cannot be used as credit towards any other product or service.

WAT042920 Revision B

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

**LC-MS System Qualification
Service Certificate**

This document covers the labor for ONE Waters® Liquid Chromatography-Mass Spectrometry System Qualification
(single or tandem quad plus one optional UV detector)

Includes: Qualification Documentation, Data System Qualification (when ordered) Organizational Binder and Qualification Overview documentation, Qualification Test Chemicals and materials as required, Qualification labor provided by a Waters Certified Compliance Specialist and Local (Zone 1) travel charges.

Pre-use Qualification Services for New Systems
Installation Qualification Documentation will be completed
Operational Qualification testing and documentation will be completed
System Performance testing will be completed and documented

Post-use Qualification Services for Installed or Existing Systems
Maintenance documentation completed for Performance Maintenance service if purchased
Operational Qualification testing and documentation will be completed
System Performance testing will be completed and documented

This certificate is for documentation purposes only and cannot be used as credit towards any other product or service.

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Waters

One-Time ACQUITY UPLC™ System Qualification Service Certificate

Part Number 741000216

Includes:

- Qualification Protocols for ACQUITY UPLC System.
- Data System Qualification (Empower or MassLynx).
- Organizational Binder and Qualification Overview documentation.
- Qualification test chemicals and materials as required.
- Qualification service labor provided by a Waters Certified Compliance Service Specialist.
- Local (Zone 1) travel charges.
- Qualification Service For New Systems
 - Installation Qualification and Operational Qualification testing and documentation will be completed.
 - System Performance testing is completed and documented for systems with ACQUITY UPLC Detector.
- Qualification Service For Installed / Existing Systems
 - A Performance Maintenance Kit must be installed (additional billed service) before Maintenance documentation can be completed.
 - Operational Qualification testing and documentation will be completed.
 - System Performance testing is completed and documented for systems with ACQUITY UPLC Detector.

*This certificate is for documentation purposes only and cannot be used as credit towards other products or services from Waters Corporation.

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The screenshot shows the 'System Strategy' dialog box with the 'Other' tab selected. The 'Other Strategy' section contains several options: 'Application Timeout (A):' set to 60 minutes, 'Disallow Use of Annotation Tools' (unchecked), and 'Permanent ID in Project' (unchecked). The 'Date Display Strategy' section has a sub-section for '(GXP/ER/ES) Time Zone Format' with 'Display Region Name' (unchecked) and 'Display Region Abbreviation' (checked). The 'Use Long Date Format (U)' checkbox is checked, and 'Display Time Zone Offset' is unchecked. The 'Example Date Format' shows '2019年04月24日 17:47:33 HKT'. The 'Advanced Database Performance' section has 'Pre-extraction Rows' set to 1000.

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Sign Off Method Name	Sign Off Method Id	Sign Off Method Version
IQ All Node summary	3201	5
IQ All Node summary	3201	5

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	结果组名	结果组日期	系统名称
1	Internal Std Sample Set	2019/5/13 18:36:01 HKT	486_sys
2	External Std Sample Set	2019/5/13 18:35:57 HKT	Alliance
3	Internal Std Sample Set	2019/5/13 18:22:11 HKT	486_sys
4	External Std Sample Set	2019/5/13 18:22:04 HKT	Alliance
5	Internal Std Sample Set	2019/5/13 17:53:46 HKT	486_sys
6	External Std Sample Set	2019/5/13 17:53:42 HKT	Alliance
7	Internal Std Sample Set	2019/5/13 17:53:23 HKT	486_sys
8	External Std Sample Set	2019/5/13 17:53:16 HKT	Alliance

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- 阅读设置和操作检定步骤：
 - SystemsQT Instrument and System Qualification Setup Procedure for the 2695/2695D Alliance and 2795 Alliance HT System with 2487/2489/2996/2998 Detectors (《采用 2487/2489/2996/2998 检测器的 2695/2695D Alliance 和 2795 Alliance HT 系统的 SystemsQT 仪器与系统检定设置步骤》)
 - Operational Qualification Procedure for the 2695/2695D/2795 Sample Heater/Cooler, Column Heater and Column Heater/Cooler (《2695/2695D/2795 样品加热器/冷却器、柱温箱和带制冷功能的柱温箱的操作检定步骤》)
 - Manual Operational Qualification Procedure for e2695 Flow Rate Accuracy Test (《e2695 流速准确度测试的手动操作检定步骤》)
- 同时设置所有温度区域
 - 每次测量间无需等待

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DETECTOR RESPONSE CHECK REPORT

ACQUISITION INFORMATION

System Name: e2695_2489 MSG
Acquired By: Melissa
Project Path: RC2_SQT_e2695_2489 final

Sample Set Id: 6422
Result Set Id: 6427

The Detector response check injection is a 0.22 mg/mL Caffeine solution injected at 7.0 μ L once the system and column have been equilibrated.

The observed height of the Caffeine peak from the detector response check injection is used to scale up the injection volume of System Readiness check and Detector Linearity & Sensitivity such that the height of the System Readiness check and Detector linearity & Sensitivity injection will be at least 1.80 AU.

This test is only for informational and diagnostic purpose. No acceptance criteria are associated with the detector response check injection.

Injection Volume (μ L) for Linear Detector Response : 13

The relationship between injection volume and peak height is approximately one-to-one. (e.g. 10% increase in injection volume results in 10% increase in peak height). Using this relationship between the injection volume and peak height for the detector response check, an injection volume as mentioned above should produce a Caffeine peak with a height of approximately 1.80 AU.

Chromatographic Peak Results

Sample Name	Injection Number	Injection Volume (μ L)	Retention Time (min)	Height (μ V)	Result Id
1 Caffeine 0.22 mg/mL	1	7.00	1.21	1104008	6437


DETECTOR RESPONSE CHECK REPORT

ACQUISITION INFORMATION

System Name: e2695_2489 MSG
Acquired By: Melissa
Project Path: RC2_SQT_e2695_2489 final

Sample Set Id: 6422
Result Set Id: 6427

Caffeine Peak for 7.0 μ L injection volume of 0.22 mg/mL solution

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
Project Path: RC2_SQT_e2695_2489 final

Processing Method: 01a System Readiness Check PM

Sample Kit Expiration Date: Friday, August 16, 2019 12:00:00 AM EDT

Sample Kit Lot Number: W16051621

Report Method: 01a System Readiness Check RM

Acquired By: Melissa

Sample Set Name: System Readiness Check SSM

Column Serial Number: 01833630013817

Sample Set Id: 6478

Result Set Id: 6484

The System Readiness Check is performed prior to the qualification to demonstrate the system and column are in an equilibrated, ready state for qualification.

The system is considered to be in a ready state only when the resulting data of Caffeine has an Area %RSD < or = 1.0, Height %RSD < or = 1.5, Retention Time %RSD < or = 1.0, Uracil peak Retention Time range is in between 0.47 to 0.95 minutes and the height of Caffeine peak is in between 1.80 AU to 2.40 AU.


Results are rounded to the same precision as displayed by the acceptance criterion prior to comparing against the specification.

Caffeine Peak Area and Height Results						
Area Spec	Area (%RSD)	Area System Status	Height Spec	Height (%RSD)	Height System Status	
1	%RSD < or = 1.0	0.0	System Ready	%RSD < or = 1.5	0.0	System Ready

Caffeine Peak RT Results			
Retention Time Spec	Retention Time (%RSD)	Retention Time System Status	
1	%RSD < or = 1.0	0.1	System Ready

Uracil Peak RT Range Results			
Retention Time Range Spec	Average Retention Time (minutes)	Retention Time System Status	
1	0.47 to 0.95 minutes	0.64	System Ready

Detector Response (Caffeine Peak) Results						
Sample Name	Detector Response Spec	Injection Volume (µL)	Average Height (AU)	Average Height (AU)	Detector Response System Status	
1	Caffeine 0.22 mg/mL	1.80 to 2.40 AU	13.00	189483.20	1.89	System Ready

SYSTEMSQT

SYSTEM READINESS CHECK REPORT

ACQUISITION INFORMATION

System Name: e2695_2489 MSG

Project Path: RC2_SQT_e2695_2489 final

Sample Set Id: 6478

Result Set Id: 6484

System Readiness Check Caffeine Peak Results

Peak Name	Val	Injection Number	Injection Volume (µL)	Retention Time (min)	Area (µV*min)	Height (µV)	Channel	Result Id
1	Caffeine	100	1	13.00	1.22	827307.40	189483.02	VC489 CHA 9495
2	Caffeine	100	2	13.00	1.22	8271603.83	1895007.37	VC489 CHA 9496
3	Caffeine	100	3	13.00	1.22	8272577.80	1895171.31	VC489 CHA 9497

System Readiness Check Uracil Peak Results

Peak Name	Val	Injection Number	Injection Volume (µL)	Retention Time (min)	Area (µV*min)	Height (µV)	Channel	Result Id
1	Uracil	100	1	13.00	0.04	981380.02	239223.91	VC489 CHA 9495
2	Uracil	100	2	13.00	0.04	980910.21	239199.28	VC489 CHA 9496
3	Uracil	100	3	13.00	0.04	981421.88	239045.71	VC489 CHA 9497

SYSTEM READINESS CHECK SUMMARY RESULTS

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QUALIFICATION SUMMARY REPORT

System Name: Qual_Training_Two
Processing Node: Waters-pcbbmm
Project Path: Qual_Training_Two_11Dec2017

Automated Qualification Tests Summary

Qualification Test Parameters	Qualification Test Specifications	Qualification Test Pass / Fail
1 System Precision Area	%RSD < or = 1.0	PASS
2 System Precision Height	%RSD < or = 1.5	PASS
3 System Precision Retention Time	%RSD < or = 1.0	PASS
4 Detector Linearity Area	R ² > or = 0.990	PASS
5 Detector Linearity Height	R ² > or = 0.990	PASS
6 Detector Area Sensitivity	%RSD < or = 4.50	PASS
7 Detector Height Sensitivity	%RSD < or = 4.50	PASS
8 Injector Linearity	R ² > or = 0.9990	PASS
9 Injector Accuracy	< or = 1.0 µL	PASS
10 Injector Carry Over	< or = 0.01 %	PASS
11 Flow Rate Linearity	R ² > or = 0.90	PASS
12 Flow Rate Accuracy	< or = 100 µL/min	PASS
13 Compositional Precision	RT %RSD < or = 3.0	Not All Peaks Found
14 Detector Noise	< or = 60 µAU	PASS
15 Detector Drift	< or = 10 mAU/hr	PASS
16 Signal to Noise	> 6000	Not All Peaks Found
17 Wavelength Accuracy 1	205nm to 207nm	PASS
18 Wavelength Accuracy 2	271nm to 275nm	PASS
19 Compositional Precision	RT %RSD < or = 3.0	PASS
20 Signal to Noise	> 6000	PASS

Manual Qualification Tests Summary

Qualification Test Parameters	Qualification Test Specifications	Qualification Test Pass / Fail
1 Sample Temperature Set Point 1	± 0.0 °C of Set-Point Temperature	PASS
2 Sample Temperature Set Point 2	± 0.0 °C of Set-Point Temperature	PASS
3 Column Temperature Set Point 1	± 0.0 °C of Set-Point Temperature	PASS
4 Column Temperature Set Point 2	± 0.0 °C of Set-Point Temperature	PASS

Waters Project Name and Version: SQT_E_2695_2489 v4.00
Date: Tuesday, January 2, 2018
Printed: 2:56:41 PM US Eastern

© 2007 - 2017 Waters Corporation. All rights reserved.

Page: 1 of 2

QUALIFICATION SUMMARY REPORT

System Name: Qual_Training_Two
Processing Node: Waters-pcbbmm
Project Path: Qual_Training_Two_11Dec2017

Automated Qualification Sample set details

Sample Set Id	Result Set Id	Sample Set Name	Sample Set Start Date
1	6152	System Precision 250UL SSM	Tuesday, December 12, 2017 3:07:42 PM EST
2	6155	Wavelength Accuracy SSM	Tuesday, December 12, 2017 3:33:36 PM EST
3	6254	Detector Lin. Sens. SSM	Tuesday, December 12, 2017 4:12:36 PM EST
4	6292	Injector Lin. Acc. SSM	Tuesday, December 12, 2017 4:27:24 PM EST
5	6302	Injector Carryover SSM	Tuesday, December 12, 2017 4:46:35 PM EST
6	6375	Flow Rate Lin. Acc. SSM	Tuesday, December 12, 2017 5:23:54 PM EST
7	6404	Compositional Prec. SSM	Tuesday, December 12, 2017 5:52:16 PM EST
8	6452	Noise and Drift SSM	Tuesday, December 12, 2017 6:46:46 PM EST
9	6472	Signal to Noise SSM	Tuesday, December 12, 2017 7:12:07 PM EST
10	6509	Wavelength Accuracy SSM	Wednesday, December 13, 2017 9:21:01 AM EST
11	6567	Compositional Prec. SSM	Wednesday, December 13, 2017 10:00:11 AM EST
12	6675	Signal to Noise SSM	Wednesday, December 13, 2017 11:30:30 AM EST

Manual Qualification Sample set details

Sample Set Id	Result Set Id	Sample Set Name	Sample Set Start Date
1	6060	SVC Temperature Acc. SSM	Tuesday, December 12, 2017 2:18:14 PM EST
2	6080	CV Temperature Acc. SSM	Tuesday, December 12, 2017 2:18:51 PM EST

Waters Project Name and Version: SQT_E_2695_2489 v4.00
Date: Tuesday, January 2, 2018
Printed: 2:56:41 PM US Eastern

© 2007 - 2017 Waters Corporation. All rights reserved.

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Sample management systems/ injectors	Injector assembly	Run instrument OQ and run system suitability
	Syringe drive assembly/syringe pump	Run instrument OQ and run system suitability
	Sample positioning drive	Verify position location and run system suitability
	Pressure transducer	Pressure test and run system suitability
	Needle/needle carriages/seal pack	Pressure test and run system suitability
	Valve seals	Pressure test and run system suitability
	Metering syringe drive assembly	Run instrument OQ

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
Self Service – License Registration

Menu | | ◀ Save as Variant... Back Exit Cancel System ▶ Execute Get Variant...

Selection-Criteria

User ID:

License Software Type:* ▼

License Status: ▼

Double Click Here To See Instructions With Images In UC:

This function will allow you to conduct the following:

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Self Service – License Registration

Menu Save Back Exit Cancel System Sort Ascending Sort descending Set filter Local file... Change layout... Select All Deselect All

Self Service – License Registration

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Self Service – License Registration

Menu Save Back Exit Cancel System Sort Ascending Sort descending Set filter Local file... Change layout... Select All Deselected All Register Serial No.

Self Service – License Registration

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Waters Licensing Activation Center

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Welcome to the Waters License Activation Center

You will need your software license serial numbers and your Software Registration file to activate your licenses or your License Deactivation file to deactivate your licenses. Need Assistance? [Contact your local office.](#)

Please select

- ☐ Breeze 2
- ☐ Empower 2
- ☒ Empower 3
- ☐ UNIFI
- ☐ NuGenesis
- ☐ Paradigm Scientific Search
- ☐ Symphony

Please select

- ☒ Workstation
- ☐ Workgroup
- ☐ Enterprise

Please select

- ☒ Activate License(s)
- ☐ Deactivate License(s)

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The screenshot shows a web browser window with the URL www.waters.com/activate/licenseintro.htm. The page header includes the Waters logo and the tagline "THE SCIENCE OF WHAT'S POSSIBLE.™". A user greeting "Welcome Chris Proehl" and "United States" is visible. The main section is titled "Please Verify/Update Your Contact Information" and includes a note: "Waters does not share your information, view our [Privacy](#) policy. Required fields are marked with an asterisks (*)." The form contains the following fields:

- * Salutation: Select Prefix (dropdown menu) with a "Clear" link.
- * First Name: text input field.
- * Last Name: text input field.
- * Email Address: text input field.
- * Company: text input field.
- * Telephone Number: text input field.
- Extension: text input field.
- * Street Address: text input field.
- * City: text input field.
- * State/Province/Region: New Jersey (dropdown menu).
- * Zip/Postal Code: 01757 (text input field).
- * Country: United States (dropdown menu).

Below the contact information fields is a section titled "Please Enter Your Order Number or Software Serial Number(s)" with the instruction: "Enter your order number to select your activations from a list of all your software serial numbers. Or, enter individually each software serial number to activate (these are displayed on the license certificate or the original software media). Your activation must include a base license if one has not already been activated." This section includes an "Order Number:" text input field and a green "Next" button.

At the bottom of the browser window, a file path "2017050411190627...lic" is visible in the address bar, and a "Show all" button is present in the bottom right corner.

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The screenshot shows a web browser window with the URL `www.waters.com/activate/licenseintro.htm`. The page header includes the Waters logo and the tagline "THE SCIENCE OF WHAT'S POSSIBLE.™". A user greeting "Welcome Chris Proehl" and "United States" is visible. The main section is titled "Please Verify/Update Your Contact Information" and includes a privacy policy notice. Below this is a form with the following fields: Salutation (dropdown menu), First Name, Last Name, Email Address, Company, Telephone Number, Extension, Street Address, City, State/Province/Region (dropdown menu showing "New Jersey"), Zip/Postal Code, and Country (dropdown menu showing "United States"). A link "Clear contact information" is next to the Salutation field. Below the contact form is a section titled "Please Enter Your Order Number or Software Serial Number(s)" with instructions and a text input field for the Order Number. A green "Next" button is at the bottom right of the form area. The browser's taskbar at the bottom shows a file named "2017050411190627...lic" and a "Show all" button.

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Welcome Chris Proehl United States

Please Verify/Update Your Contact Information

Waters does not share your information, view our [Privacy](#) policy. Required fields are marked with an asterisks (*).

* Salutation: Select Prefix [Clear contact information](#)

* First Name:

* Last Name:

* Email Address:

* Company:

* Telephone Number:

Extension:

* Street Address:

* City:

* State/Province/Region: New Jersey

* Zip/Postal Code: 01757

* Country: United States

Please Enter Your Order Number or Software Serial Number(s)

Enter your order number to select your activations from a list of all your software serial numbers. Or, enter individually each software serial number to activate (these are displayed on the license certificate or the original software media). Your activation must include a base license if one has not already been activated.

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Select Software Serial Number(s) to Activate

Base License Must Be Activated First

If your base license has not been activated you need to (1) select it from this order along with any other licenses you would like to activate, or (2) if it does not appear on this order, go back and activate it first, then proceed to activating the licenses on this order.

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Please Verify/Update Your Contact Information

Waters does not share your information, view our [Privacy](#) policy. Required fields are marked with an asterisks (*).

The following error(s) have occurred.

Need Assistance? [Contact your local office.](#)

- EMP011001A00509 - Invalid Serial Number. Please verify you have entered the correct serial number.
- EMP011003A00402 - Serial number not valid for product indicated.
- W6DNL8047M - Serial Number not valid for Company indicated. Please verify the spelling of your company name.

* Salutation:	Mrs.	Clear contact information
*First Name:	Jodie	
*Last Name:	Bouthillette	
*Email Address:	Jodie_L_Bouthillette@waters.com	
*Company:	Waters Corporation	
*Telephone Number:	508 482 8678	

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

The screenshot shows the 'Waters License Activation Center' web page. At the top is the Waters logo and tagline. Below is a welcome message and instructions to provide software license serial numbers and a Software Registration file. There are three sections for selection, each with radio buttons. The first section, 'Please select', lists software products: Breeze 2, Empower 2 (selected), Empower 3, UNIFI, NuGenesis, and Paradigm Scientific Search. The second section, 'Please select', lists configurations: Workstation (selected), Workgroup, and Enterprise. The third section, 'Please select', lists actions: Activate License(s) and Deactivate License(s) (selected). A green 'Next' button is at the bottom.

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Please select

- ☐ Activate License(s)
- ☒ Deactivate License(s)

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Load and Deactivate Your License Deactivation File

Select the License Deactivation file you generated from the Waters Licensing Wizard. You will need your License Deactivation file to deactivate your licenses. Need Assistance? [Contact your local office.](#)

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