

Customer:

Dr. Qingsong Lin
National University of Singapore
Department of Biological Sciences
14 Science Drive 4
Singapore 117543

Email : dbslinqs@nus.edu.sg
Tel : 65 6874 1682
Fax : 65 6872 2013

Quotation for MS

Quotation No. 22147814 - Expiration Date : 30.01.2020

General information

111126/10263420

Provider:

Waters Pacific Pte Ltd
Julius Hiu
1 Science Park Road
#02-1/06 The Capricorn
Singapore Science Park II
Singapore 117528

Tel : 65 83219687
Email : Julius_Hiu@waters.com

Quotation
Quotation for MS

Price Details in SGD

Please reference this quotation when purchase order is issued

Itm	Part No.	Description	Qty	Net Price
		Xevo G2-XS System	1	
		<i>With the following configuration:</i>		
1	176003480	Xevo G2-XS QToF 4k	1	
2	176002846	Vacuum Pump	1	
3	176002662	nanoLockSpray source w/universal sprayer	1	
4	176002507	NanoFlow Capability Kit	1	
5	186007529	LOW FLOW ESI PROBE	1	
6	176002525	MassLynx Performance Workstation	1	
7	668000273	MONITOR, Lenovo ThinkVision Flat Panel	1	
8	176002513	MS Ref Stds: TOF G2-S MS Install	1	
9	186009340	Nitrogen Generator, Genius XE35, 230V	1	
10	176016013	ACQ UPLC M-Class:uBSM, uSM-FL, TVM, ASM	1	
11	186007472	ACQ M-Class HSS T3 1.8µm 300µmx150mm	1	
12	205001681	Kit, ACQ M-Class 75/300um Scale ZenFit	1	
13	186002325	MassPrep Enolase Digestion Standard	1	
14	205001681	Kit, ACQ M-Class 75/300um Scale ZenFit	1	
15	741000364	XEVO G2 SYSTEM INSTALL CERTIFICATE	1	
16	750000840	5-day On-Site MS App Training w/eLearn	1	
17	176003951	Nano/Micro LC-MS Solvent Install Kit	1	
18	176001086	TargetLynx	1	
19	176004040	Progenesis QI & QI for proteomics	1	
20	176004029	Progen Proteomics High Throughput Wrkstn	1	
21	668000273	MONITOR, Lenovo ThinkVision Flat Panel	1	
22	176003497	EZ Info	1	
		Waters Xevo TQ-XS System	1	
		<i>With the following configuration:</i>		
23	176003945	Waters Xevo TQ-XS	1	
24	176002526	MassLynx Workstation with TL	1	
25	668000273	MONITOR, Lenovo ThinkVision Flat Panel	1	
26	176003946	Dry Vacuum EBARA pump	1	
27	176003948	MS Ref StdsXEVO TQ-XS	1	
28	186007529	LOW FLOW ESI PROBE	1	
29	176016013	ACQ UPLC M-Class:uBSM, uSM-FL, TVM, ASM	1	
30	186007472	ACQ M-Class HSS T3 1.8µm 300µmx150mm	1	
31	205001681	Kit, ACQ M-Class 75/300um Scale ZenFit	1	
32	186002325	MassPrep Enolase Digestion Standard	1	
33	205001681	Kit, ACQ M-Class 75/300um Scale ZenFit	1	
34	176001646	Nano-flow ESI	1	
35	741000358	XEVO Tandem Quad System Install Cert.	1	

Quotation
Quotation for MS

Price Details in SGD

Please reference this quotation when purchase order is issued

Itm	Part No.	Description	Qty	Net Price
36	750000530	On-Site MS/App Training -ROW (3 Days)	1	
37	176003951	Nano/Micro LC-MS Solvent Install Kit	1	
38	740005170	TAW Xevo G2-XS QToF (1PM)	1	
39	740005074	TAW ACQUITY M-Class µBSM (1PM)	2	
40	740005076	TAW ACQUITY M-Class µSM-FL (1PM)	2	
41	740005080	TAW ACQUITY M-Class TVM (0PM)	2	
42	740005078	TAW ACQUITY M-Class ASM (1PM)	2	
43	740005659	TAW Xevo TQ-XS (1PM)	1	
44	740005082	MAINT - Progenesis QI	1	
45	740005084	MAINT -Progenesis QI for Proteomics	1	
46	740002549	Maint: MassLynx S/W 2nd Yr	1	

Total excluding GST			1,089,320.85
GST	7.00 %		76,252.46
Total			1,165,573.31

Julius Hiu
Sales Manager, Singapore

Mathieu Laouenan
General Manager - SEA & Australasia

Terms of Payment

30 Days - document date

Shipping Information

DDP Singapore

Quotation

Quotation for MS

Optional Products

Please reference this quotation when purchase order is issued

Itm	Part No.	Description	Qty	Net Price
47	176003403	ionKey Source w/Infill Assy - Xevo TQ-S	1	32,279.40
48	176003949	APCI Probe Option	1	2,798.55
49	205016040	ACQUITY UPLC M-Class Cart	1	2,255.40

Quotation

Quotation for MS

Detail Product Descriptions

Item	Part No.	Description
------	----------	-------------

1	176003480	Xevo G2-XS QToF 4k
---	-----------	---------------------------

The Waters Xevo® G2-XS QToF is a high performance benchtop QToF mass spectrometer incorporating the XS Collision Cell for enhanced sensitivity and selectivity, resulting in the highest quality most comprehensive qualitative information together with the ultimate in quantitative performance.

For the scientist who needs to identify, quantify and confirm the broadest range of compounds in the most complex and challenging samples, the Xevo® G2-XS QToF offers maximum robustness with no compromise in performance, class-leading real-world quantitative sensitivity, flexibility to adapt to changing needs, and is accessible to experts and non-experts alike. STEPWAVE" ION OPTICS

StepWave" ion transfer optics dramatically increases the efficiency of ion transfer from the ion source to the MS analyser at the same time as efficiently eliminating undesirable neutral contaminants. This technology delivers unprecedented levels of robust, durable sensitivity.

2. INTELLISTART" FLUIDICS

The instrument is equipped with an on-board infusion system capable of delivering reference solutions from 3 built-in vial locations. The on-board fluidics system is controlled by the IntelliStart" software to provide automated instrument setup and mass calibration. The reference solutions are delivered via switching valves for either direct or combined (into an LC flow) infusion into the ion source. The valves can in addition be programmed from the software to function as an LC flow divert. If required, the fluidics can be controlled manually via the System Console.

3. QUADRUPOLE ANALYSER

The instrument is equipped with a high performance quadrupole mass analyser. Pre-filters are fitted to the mass analyser to maximise resolution and transmission. The pre-filters also eliminate the need for cleaning of the quadrupole mass analyser. All lens and analyser voltages are digitally controlled. Analyser parameters may be programmed with respect to mass for optimal performance.

4. COLLISION CELL

The instrument is equipped with an XS Collision Cell where the collision energy is under programmable software control. The collision energy and collision gas pressure used during data acquisition are automatically recorded and appended to the relevant data file.

5. TOF ANALYSER

The instrument is equipped with a high performance oaToF mass analyser with a mass range up to m/z 100,000 and a resolving power of 40,000 FWHM.

6. VACUUM SYSTEM

The instrument is equipped with a differentially pumped, automated vacuum system comprising air-cooled turbomolecular pumps and one backing pump (either one rotary pump or one oil free pump). Vacuum read backs and system vent/pump cycles are digitally monitored and controlled, to provide total software control and ensure fail-safe operation in the event of power failure.

The backing option must be ordered separately (See Backing Pump Options for part numbers and descriptions)

LockSpray " API interface#Dual orthogonal interface for robust LC/MS is included.

The standard system does not include the following items, which must be specified separately:

nanoLockSpray" interface for nanoflow chromatography

Additional inlet probes and ion source options (detailed below).

Acquisition PC data system and monitor

Additional Workstation data system terminals.

Printers.

Vacuum backing pump options (rotary pump or oil-free pump)

Quotation

Quotation for MS

Detail Product Descriptions

Itm	Part No.	Description
------------	-----------------	--------------------

HPLC systems or other inlet options.

2	176002846	Vacuum Pump
----------	------------------	--------------------

3	176002662	nanoLockSpray source w/universal sprayer
----------	------------------	---

nanoLockSpray source w/universal sprayer

This is a dual nanospray inlet, which is able to supply a reference compound separately from the analyte spray; selection between the sprays is achieved via an indexed rotating baffle located in front of the sample cone. This avoids problems associated with post column addition of a reference compound such as ionisation suppression, mass interference and solvent gradient effects.

The nanoLockSpray source comes complete with the Universal NanoFlow sprayer and camera for visualisation of nanoLockSpray on the acquisition PC.

4	176002507	NanoFlow Capability Kit
----------	------------------	--------------------------------

NanoFlow Capability Kit

Idex flow controller for all nanoFlow applications.

Note: This kit is required for use with the nanoLockspray and Trizaic source options.

5	186007529	LOW FLOW ESI PROBE
----------	------------------	---------------------------

6	176002525	MassLynx Performance Workstation
----------	------------------	---

MassLynx Application Software License and License for OpenLynx"

MassLynx OpenLynx

Summary reports for groups of samples. OpenLynx methods may be configured to include:

1) Analytes found/not found. 2) Purity, 3) Chromatograms, 4) Spectra. OpenLynx reports are portable and may be viewed (read-only) with optional Browser utility without MassLynx software. Mass Spectrometer Control Software

Contains the performance MS workstation as defined in the MassLynx Variant configurator 176706000 and Product Ordering Guide

8	176002513	MS Ref Stds: TOF G2-S MS Install
----------	------------------	---

9	186009340	Nitrogen Generator, Genius XE35, 230V
----------	------------------	--

Quotation

Quotation for MS

Detail Product Descriptions

Item	Part No.	Description
------	----------	-------------

10	176016013	ACQ UPLC M-Class:uBSM, uSM-FL, TVM, ASM
-----------	------------------	--

ACQ UPLC M-Class:uBSM, uSM-FL, TVM, ASM

The ACQUITY UPLC M-Class Core System has solvent, sample, auxiliary, and trap valve managers. The system is configured for micro-scale binary solvent delivery (¼BSM) and sample introduction with a fixed loop sample manager (¼SM-FL).

The following items are included as part of this system configuration:

- Micro Binary Solvent Manager (¼BSM, 186016003)
- Micro Sample Manager # Fixed Loop (¼SM-FL, 186016012)
- Auxiliary Solvent Manager (ASM, 186016022)
- Startup Kit, ACQUITY UPLC M-Class System, with Solvent Tray Module (200000405)
- Trap Valve Manager (TVM,186016016)

Supported Configurations:

- direct injection
- single pump trapping (trap/elute)
- dual pump trapping (trap/elute) with or without ASM The ASM

The ASM delivers two (2) separate isocratic eluents; one for nanolockspray addition and another that can be used as an additional trapping eluent or other isocratic application. The ¼BSM can be used for single pump trapping or the ASM can be used to deliver an alternate trapping eluent.

11	186007472	ACQ M-Class HSS T3 1.8µm 300µmx150mm
-----------	------------------	---

ACQUITY UPLC M-Class HSS T3 Column, 100Å, 1.8 µm, 300 µm X 150 mm, 1/pk

12	205001681	Kit, ACQ M-Class 75/300um Scale ZenFit
-----------	------------------	---

13	186002325	MassPrep Enolase Digestion Standard
-----------	------------------	--

MassPREP Enolase Digestion Standard

15	741000364	XEVO G2 SYSTEM INSTALL CERTIFICATE
-----------	------------------	---

Xevo Q-TOF System Installation

Includes:

- System Set up and Specification Testing
- Product Familiarization Training
- 1 Year Manufacturers Warranty
- Insight I-Assist Remote Connection

Quotation

Quotation for MS

Detail Product Descriptions

Item Part No.	Description
16 750000840	5-day On-Site MS App Training w/eLearn
17 176003951	Nano/Micro LC-MS Solvent Install Kit
	Nano/Micro LC-MS Solvent Install Kit The Nano/Micro LC-MS Solvent Install Kit is intended to assist in effectively installing new nano or micro flow LC-MS systems by minimizing sources of possible contamination. This kit includes the LC- MS Grade Solvent Kit (186008715) as well as the Waters Certified Container Low Volume Kit (186007278). Additional kit components also include a formic acid additive, aqueous ammonia additive and nitrile gloves. Note: Because this is a hazardous material and GHS requirements vary from country to country, Waters Corporation is only able to ship this item to specific countries worldwide. If the destination for this order (i.e. the ship to country) is changed to a non-permitted country at any point after this quote was generated, this item will be excluded from the final order.
18 176001086	TargetLynx
	The TargetLynx Application Manager is designed for target quantification. It includes all the functionality available in QuanLynx plus comprehensive quality control (QC) features. The identification of target analytes is assured based on user-defined ion ratios. Quality limits are user definable; data failing to meet the defined criteria are automatically flagged. In addition TargetLynx highlights significant results where either the 'surveillance' reporting level or action level is exceeded. QC parameters available include: - Up to six ion ratios with automated batch to batch up-date. - Retention time tolerance with automated batch to batch up-date. - Signal to noise ratio. % recovery calculation. - Calculation of variation in response factor. TargetLynx fulfills analyte confirmation requirements in food analysis when analyzing for regulated pesticide or antibiotic residues. Additionally, TargetLynx is designed for confirmation of dioxin assays. When installed on either a Waters 3100 Mass Detector, SQ Detector, TQ Detector, Xevo TQ MS or Xevo TQ-S MS, TargetLynx also includes Quanpedia, QC Monitor and TrendPlot, providing you with compound and method database functionality, pro-active quantitative data quality checking as well as intra- and inter-batch data and performance checking.
19 176004040	Progenesis QI & QI for proteomics
20 176004029	Progen Proteomics High Throughput Wrkstn
21 668000273	MONITOR, Lenovo ThinkVision Flat Panel
	Monitor, Lenovo Thinkvision Flat Panel 22 inch monitor

Quotation

Quotation for MS

Detail Product Descriptions

Item	Part No.	Description
------	----------	-------------

22	176003497	EZ Info
----	-----------	---------

23	176003945	Waters Xevo TQ-XS
----	-----------	-------------------

Xevo® TQ-XS is an advanced benchtop tandem quadrupole mass spectrometer which expands the scope of ultimate sensitivity analysis. It features StepWave XS™ ion transfer optics delivering enhanced sensitivity, robustness, and reliability; passively removing gas load and undesirable neutral contaminants while actively transferring ions into the mass analyzer. Sensitivity benefits are more accessible as the Xevo TQ-XS also features a tool-free probe design which allows users much easier maintenance, optimization and improved reproducibility between operators. Method transfer onto the Xevo TQ-XS is made simple by the Extended Dynamic Range (XDR™) detector which allows six orders of linear dynamic range. Wider compound coverage without changing ionization technique is delivered by the revolutionary UniSpray™ source option. In addition to these benefits, users can combine quantitative MRM with qualitative MS data acquisition through RADAR™ and confirm low levels of compounds using ScanWave™ enhanced Product Ion Confirmation scanning. The MassLynx® software supporting this instrument ensures that automated system optimization, status monitoring and results review and reporting is a smooth process using IntelliStart™ and TargetLynx XS™. When combined these advances allow new levels of reproducible and accurate quantitation of analytes at the lowest concentrations.

The following items are included as part of the standard system:

- Z SPRAY™ API interface#Dual orthogonal interface for robust LC/MS Tool-free electrospray (ESI) inlet probe#for efficient ionisation of a wide range of compounds
- ESCI™ ionisation capability#rapid switching source for both ESI and APCI in the same run
- IntelliStart™ fluidics#Automated tuning, calibration and method development
- TargetLynx XS™#Application manager (requires license as provided with PC, below)
- OpenLynx™#Application manager (requires license as provided with PC, below)

The standard system does not include the following items, which must be specified separately:

- Additional inlet probes and ion source options (detailed below).
- Acquisition PC data system and monitor
- Additional Workstation data system terminals.
- Printers.
- Additional MassLynx™ options (detailed below).
- Vacuum backing pump options (rotary or oil-free combinations)
- HPLC systems or other inlet options.

1. Z SPRAY™ API INTERFACE

This instrument is equipped with an atmospheric pressure ionisation (API) LC interface. The source and spraying elements are visible through a transparent window in the enclosure and are easily accessible via a quick-release mechanism. The source elements may be wiped clean in situ or removed for cleaning without the need for tools and without breaking vacuum. The nebulized spray is orientated orthogonally and positioned off axis for maximum source longevity and analyser protection against 'dirty' samples. The source also includes facilities for de-clustering ions formed at atmospheric pressure. Positive and negative capability is included. Positive ion, negative ion and ESCI™ capability is available as standard (allowing rapid switching between ESI and APCI, positive and negative in the same run). All source voltages and gases are under data system

Quotation

Quotation for MS

Detail Product Descriptions

Item	Part No.	Description
------	----------	-------------

control. Both ESI and APCI probes are now tool-free allowing easy routine maintenance and improved user-to-user reproducibility.

2. INTELLISTART™ FLUIDICS

The instrument is equipped with an on-board infusion system capable of delivering reference solutions from 3 built-in vial locations. The on-board fluidics system is controlled by the IntelliStart software to provide automated instrument setup, mass calibration and method development. The reference solutions are delivered via switching valves for either direct or combined (into an LC flow) infusion into the API source. The valves can in addition be programmed from the software to function as an LC flow divert. If required, the fluidics can be controlled manually via the system Console.

3. TANDEM QUADRUPOLE ANALYSER

The instrument is equipped with two high performance quadrupole mass analysers with inter-element beam focusing and a mass range of 2-2000 amu. Pre-filters are fitted to the mass analysers to maximise resolution and transmission. The pre-filters also eliminate the need for cleaning of the quadrupole mass analysers. All lens and analyser voltages are digitally controlled. Analyser parameters may be programmed with respect to mass for optimal performance. Analyser parameters used for data acquisition are automatically recorded and appended to the relevant data file.

4. StepWave XS ION TRANSFER OPTICS

This instrument is equipped with patented (US patent - US9035241), off-axis StepWave™ XS ion transfer optics which provide improved sensitivity that is reproducible for multiple injections and can deliver reliable quantitation at the very lowest levels. StepWave XS ensures neutrals and gas load are passively removed for enhanced transmission, with the ions actively transferred into the mass analyser, improving sensitivity and robustness. The first stage is modified with a series of horizontal plates that allows a more controlled extraction of ions from the gas flow so that ions collide with gas molecules with less energy. A segmented quadrupole second stage focuses the ion beam to provide more efficient transmission of ions from the source to the detector.

5. COLLISION CELL

The instrument is equipped with a unique design of collision cell that can be operated in two different modes (T-Wave or ScanWave), depending upon need. In both modes of operation the collision energy and gas pressure are under programmable software control. The collision energy and collision gas pressure used during data acquisition are automatically recorded and appended to the relevant data file.

T-Wave™:

The cell can be operated as a high efficiency travelling wave (T-Wave) device for collision induced dissociation. The travelling wave enables rapid cell clearance and refill for fast MRM transition switching while maintaining optimum signal to noise.

ScanWave™:

The cell can be operated in a ScanWave mode. The cell can accumulate ions and then release them, according to their mass, in synchrony with the second quadrupole mass analyser. This mode of operation results in a significant increase in the signal intensity of full scan spectra (Product ion).

RADAR™:

An information-rich acquisition approach that allows you to collect highly specific quantitative MRM data for target compounds while providing additional spectral data to help visualize all other components in the sample.

6. VACUUM SYSTEM

The instrument is equipped with a differentially pumped, automated vacuum system comprising: three air-cooled turbomolecular pumps and a single oil-free vacuum backing pump. Vacuum read backs and system vent/pump cycles are digitally monitored and controlled, to provide total software control and ensure fail-safe operation in the event of power failure. The backing option must be ordered separately (See Backing Pump Option for part number and description)

7. DETECTOR

The instrument is equipped with an Xtended Dynamic Range™ (XDR) low noise dynolite photomultiplier detector. The detector is positioned after the second analyzer. A High Voltage conversion dynode and phosphor are positioned at 90° off-axis to the analyser for the elimination of neutral noise. The detector features novel, integral focusing optics, which provides a detection

Quotation
Quotation for MS

Detail Product Descriptions

Itm Part No. Description

efficiency approaching 100% for single ions. The photomultiplier is enclosed in its own vacuum envelope for long life. The detector operates in both positive and negative ion mode, which can be switched rapidly under software digital control.

7. MASSLYNX SOFTWARE / MS Workstation

The MS Workstation and MassLynx 4.2 License for the application software for instrument control, data acquisition and processing must be ordered separately (See MS Workstation Variant Configurator section for part numbers and descriptions) Post acquisition processing and general data manipulation can be carried out by an additional computer workstation and software installation (See the MassLynx Process-only Workstations section of the MassLynx Variant configurator 176706000).

24 176002526 MassLynx Workstation with TL

MassLynx" Workstation with License

MassLynx Application Software License includes key discs for TargetLynx" and OpenLynx"

Contains the standard MS workstation as defined in the MassLynx Variant configurator 176706000 and Product Ordering Guide

26 176003946 Dry Vacuum EBARA pump

The Ebara EV-SA30 incorporates totally clean multi-stage Roots technology with no oil or tip seal dust contamination risk. No annual service is required and it has 80% lower power consumption and 50% smaller footprint compared to two rotary vane pumps giving a similar combined pumping speed.

27 176003948 MS Ref StdsXEVO TQ-XS

34 176001646 Nano-flow ESI

Nano-flow ESI

An interface enabling nano-flow (<1µl/min) electrospray.

Nano-flow ESI comprises:

- Universal nano-flow Sprayer.
- Built in camera.
- Adjuster assembly.

Note: It is recommended that a harvard syringe pump is ordered with this source option to help with instrument tuning. VB-HSP11 Harvard Syringe Pump

35 741000358 XEVO Tandem Quad System Install Cert.

Xevo MS System Installation

Includes:

- System Set up and Specification Testing
- Product Familiarization Training

Quotation

Quotation for MS

Detail Product Descriptions

Itm	Part No.	Description
-----	----------	-------------

- 1 Year Manufacturers Warranty
- Insight I-Assist Remote Connection

36	750000530	On-Site MS/App Training -ROW (3 Days)
-----------	------------------	--

On-Site MS/Application Training - ROW (3 Days)

38	740005170	TAW Xevo G2-XS QToF (1PM)
-----------	------------------	----------------------------------

39	740005074	TAW ACQUITY M-Class μBSM (1PM)
-----------	------------------	--

40	740005076	TAW ACQUITY M-Class μSM-FL (1PM)
-----------	------------------	--

41	740005080	TAW ACQUITY M-Class TVM (0PM)
-----------	------------------	--------------------------------------

42	740005078	TAW ACQUITY M-Class ASM (1PM)
-----------	------------------	--------------------------------------

44	740005082	MAINT - Progenesis QI
-----------	------------------	------------------------------

45	740005084	MAINT -Progenesis QI for Proteomics
-----------	------------------	--

46	740002549	Maint: MassLynx S/W 2nd Yr
-----------	------------------	-----------------------------------

Masslynx Software 2 Year Plan

MassLynx Software Maintenance Contract - Valid for 2 Years - Includes coverage according to the Waters Informatics Maintenance Plan.

47	176003403	ionKey Source w/Infill Assy - Xevo TQ-S
-----------	------------------	--

ionKey Source w/Infill Assy - Xevo TQ-S

The ionKey Source with infill assembly contains the ionKey Tandem Quad source, the info set, a cable management bracket which aligns the cables and tubing from the UPLC to the source and a Xevo TQ S panel which hides the bracket and cables. The ionKey tandem quad source is designed to accept the iKey and precisely position it into the MS. The source has an external #Stop Flow# button that is used to stop the flow and decompress the iKey before removing it from the source. The source is sealed from environmental interference's decreasing background contamination.

Quotation

Quotation for MS

Detail Product Descriptions

Itm	Part No.	Description
------------	-----------------	--------------------

48	176003949	APCI Probe Option
-----------	------------------	--------------------------

Dedicated tool-free probe option for Xevo TQ-XS atmospheric pressure chemical ionisation (APCI). The interface can be used at up to 2 ml/min without the need for flow splitting.

49	205016040	ACQUITY UPLC M-Class Cart nanoACQUITY UPLC Flex Cart
-----------	------------------	--

Quotation

Quotation for MS

General Sales Terms and Conditions

1. General
Waters Pacific Pte Ltd ("Seller") hereby agrees to sell the products to buy upon the terms and conditions set forth herewith. Seller shall not be bound by any standard or printed terms furnished by the Buyer in any of its documents, unless the Buyer specifically states in writing separately from such terms that it intends such terms to apply and the Seller acknowledges such notification in writing.
2. Price
 - a. Unless otherwise agreed at the time of sale, the price stated is ex-works and exclusive of Goods & Services Tax which shall be chargeable in addition upon supply of goods.
 - b. All quotations and estimates used by Seller are, unless otherwise stated, based on current costs and are subject to amendment on or after acceptance of order by Seller to meet any rise in such costs.
 - c. Any variation to prices quoted as a result of government taxes and levies will be for the Buyer's accounts.
 - d. Written quotations automatically expire 30 calendar days from the date issued and are subject to termination by notice within that period. All orders based on written quotations are subject to Seller's acceptance at its office in 1 Science Park Road #02-01 The Capricorn, Singapore Science Park II, Singapore 117528.
3. Delivery
 - a. Stock items are processed for prompt delivery. Whilst there is no minimum order value there is a minimum charge in respect of packing, handling and delivery. Details of these charges are available on request.
 - b. Notwithstanding that the Seller and the Buyer may agree to deliver the goods at some specified place the delivery of goods to a carrier for the purpose of transmission to the Buyer is deemed to be a delivery of the goods to the Buyer and risk in the goods shall pass accordingly at the moment of delivery to the carrier.
 - c. Inspection and testing and/or installation of the goods is not provided by Seller unless otherwise agreed at the time of sale.
 - d. Seller retains ownership of the goods the property in which shall not pass to Buyer and Buyer shall keep any goods delivered to it as bailee for and on behalf of Seller until Seller has received payment of the price of all the goods (whether or not the goods are delivered in instalments and some have been paid for by Buyer pursuant to these terms) and until such time the Buyer shall store the goods separately or in some other way ensure that the goods are readily identifiable as the property of Seller and at Seller's request either shall deliver up to the goods to Seller or shall permit Seller to enter Buyer's business or other premises to repossess the goods and subsequently re-sell them.
 - e. The provisions of this clause 3(d) shall survive the determination of the contract for whatever reason.
4. Payment
 - a. Payment terms are strictly Net Cash 30 days after delivery of the goods, provided Buyer's credit has been approved or specified in the quote. If the Buyer's credit had not been established with Seller, terms are payment in advance.
 - b. All payments made by Buyer shall be made without any deduction of deferment in respect of any disputes or claims whatsoever and/or in respect of any taxes imposed by or under the authority of any government or public authority.
 - c. Where Seller does not receive payment of any of its invoices within 30 days of the date thereof interest shall thereafter accrue on the sum due and owing to Seller at the rate of 2% per month calculated on a daily basis without prejudice to Seller's right to receive payment within such 30 day period.
 - d. Time for payment shall be of the essence.
5. Confidentiality/Intellectual Property
 - a. Any order received by Seller will be treated as confidential and will not be disclosed or publicised to any third party for any reason without the prior written consent of Buyer. Neither shall Seller use Buyer's name or other identify for advertising or publicity purposes without Buyer's prior written consent.
 - b. The sale of the goods by Seller to Buyer does not confer any right or license upon buyer to use, exploit or otherwise utilise any intellectual property rights subsisting in or relating to the goods of which Seller is a proprietor or to which Seller is otherwise entitled.
6. Specifications
The characteristics and specifications of Seller's products are documented in Seller's technical literature and are approximately only Seller reserves the right to change product specifications without prior notice.
7. Cancellation
An order accepted by Seller may be cancelled only before delivery of the goods and only with the consent of Seller upon terms that Buyer will indemnify Seller against any loss incurred.
8. Return
In no case are products to be returned without first obtaining the consent of the Seller. Only unused products as currently manufactured which have been invoiced by the Seller within 30 days will be considered for return. Products accepted for credit are subject to a 20 per cent service charge plus all transportation charges. Products built to order or not of the original manufacturer of the Seller are not subject to return for credit under any circumstances. Products must be securely packed to reach the Seller without damaged. No returns shall be authorised after 90 days following shipment to buyer.
9. Warranty
The Products and/or services shall be covered by the applicable Waters' standard warranty, a copy of which is supplied with the Products and/or services or upon request. NO OTHER WARRANTY, WHETHER EXPRESS OR IMPLIED, IS MADE WITH RESPECT TO THE PRODUCTS AND/OR SERVICES. WATERS EXPRESSLY EXCLUDES THE IMPLIED WARRANTIES OF MERCHANTABILITY AND OF FITNESS FOR A PARTICULAR PURPOSE. Any model or sample furnished to the Buyer is merely illustrative of the general types and quality of goods and does not represent that the Products will conform to Waters applicable standard warranty. WATERS SHALL NOT BE LIABLE FOR CONSEQUENTIAL, INCIDENTAL, SPECIAL OR ANY OTHER INDIRECT DAMAGES RESULTING FROM ECONOMIC LOSS OR PROPERTY DAMAGE SUSTAINED BY BUYER FROM THE USE OF ITS PRODUCTS OR SERVICES.
10. Termination
 - a. It prior to the delivery of the goods (or any instalment thereof) Seller has reason to believe that Buyer is unable to pay it's debts or if Buyer enters into liquidation (other than for the purpose of effecting a reconstruction or amalgamation) whether compulsory or voluntary or compounds with or conveys a meeting of its creditors or has a receiver appointed over all or any part of its assets or takes or suffers any similar action in consequence of a debt or ceases for any reason to carry on business, Seller may terminate the contract by notice in writing to Buyer.
 - b. Termination of the contract personate to Clause 10(a) shall not effect the accrued rights of the parties arising in any way out of such contract as of the date of termination.
11. Variation
Neither Seller or Buyer shall be bound by any variation, waiver of or an addition to, those terms and conditions except as agreed by both parties in writing.
12. Time
Seller will endeavour to dispatch the products on the promised delivery date but does not guarantee to do so.
13. Force Majeure
Waters shall have no liability for failure to perform or delay in performance, in the delivery of any and all equipment manufactured or sold by Waters including instruments, supplies, components, systems, chemistry, accessories, replacement spare parts, or any and all services provided by Waters, caused by circumstances beyond its reasonable control including, but not limited to, acts of God, acts of nature, floods, fire, explosions, war or military mobilisation, United States governmental action or inaction, request of governmental authority, delays of any kind in transportation or inability to obtain material or equipment, acts of other governments, strikes or labour disturbances.
14. Law
The contract shall be deemed to have been made in Singapore and the Parties to the contract hereby submit to the jurisdiction of Singapore courts. Singapore law shall be the proper law of the contract.