

IQ-02 TEST SCRIP DATA SHEET:

AMBIENT TEMPERATURE: 21 °C to 22 °C

Verified using

THERMOHYGROGRAPH

Serial Number : 5650

EQID : 0418EQID : 0418

Calibration Due Date: 07-Jan-2021


flex
SKW-MFG4018-00
Temperature and Humidity Monitoring Record

Rev: A
Form A

(1) Month / Year : Mar 20 (4) Equipment Used : Thermohygrograph (7) EQ ID : 0418
(2) Floor Level : 5th Floor (5) Serial Number : 5650 (8) Temperature Range : 18°C to 30°C
(3) Area : Production (6) Calibration due date : 07 Jan 21 (9) Humidity Range : 25%RH to 85%RH

Date Check	(10) Time Check	Measurement				(15) Within Spec. (Check since the last recording)	(16) Out of Control Record				(21) Name / Signature	(22) NCMR
		(11) Min. Temp. (°C)	(12) Max. Temp. (°C)	(13) Min. Humidity (%RH)	(14) Max. Humidity (%RH)		(17) Temp. (°C)	(18) Humidity (%RH)	(19) Date Out	(20) Time Out		
1	X	X	X	X	X	Yes / No	X	X	X	X	Yep him wah g	X
2	0714	21	22	51	55	(Yes) No	NA	NA	NA	NA	Yep him wah g	NA
3	0734	21	22	53	54	(Yes) No	NA	NA	NA	NA	Yep him wah g	NA
4						Yes / No						
5						Yes / No						
6						Yes / No						

RELATIVE HUMIDITY: 51 % to 55 %
Verified using
THERMOHYGROGRAPH
Serial Number : 5650
EQID : 0418EQID : 0418
Calibration Due Date: 07-Jan-2021



<u>Calibration</u>	
EQID	: 0418
S/No	: 5650
Cal Date	: 07-JAN-20
Due Date	: 07-JAN-21
Remark	: NA

ELECTRICAL SUPPLY :



Verified using :

FLUKE 110 TRUE RMS MULTIMETER

Serial Number: 85730619

EQID: 0245

Calibration Due Date: 16-Dec-2020



Minimum 4 dedicated outlets

ELECTRICAL OUTLET

TE82526 Pump
Tray Tester (AC Line)

LCD Monitor (AC
Line_

Computer Work
Station (AC Line



**AND HR 200
Electronics Balance
(AC Line)**



NETWORK AVAILABLE : Minimum of 1 dedicated network

Blue Sticker : Dedicated network for Pump Tray Area



WS-9 is a dedicated network

WORK SPACE: The system has sufficient room for the test station, UUT and service



TEST STATION BENCH : The test station bench is large and strong enough to support the station



DOCUMENT TITLE:	TEST SPECIFICATION AND PROCEDURE, ALLIANCE PUMPTRAY
SYSTEM/ASSEMBLY NAME:	ALLIANCE PUMPTRAY
PART NUMBER AFFECTED	279000419, 279000524, 289000526, 700000713, 700001030, WAT270463, WAT271078
TEST EQUIPMENT NAME:	ALLIANCE PUMPTRAY TESTER
TEST EQUIPMENT NUMBER:	TE82526
TEST SOFTWARE:	*1.0*

AUTHOR(S): MIKE COTTER

DATE: April 30, 1998

REVISED BY: MIKE COTTER

DATE: May 05, 2016

Document No.: TE82526 2-JQR
Iteration No.: 00
Attachment: 7.13 Page: 1 of 2

IQ-04: Test Software and Instrument Firmware Version



IQ-05 Test Script Data Sheet:

**Station does not have Pinch Points/ sharp edges
Glass Bottle is secured by plastic container**



Glass Bottles
secured inside the
container