

Thermo Scientific

Software Installation Qualification Report

Product: Thermo Exactive Series
Version: 2.12.0.3134
Software Installation Date/Time: 04/07/2022 11:56:49
Report Creation Date/Time: 04/07/2022 11:57:58

File and Folder Installation Actions:

Action	Attempts	Pass Counts
Install File	80	80
Uninstall File	0	0
Remove File	0	0
Retain File	0	0
Create Folder	0	0
Remove Folder	0	0
Verify File	63	63

Registry Installation Actions:

Action	Attempts	Pass Counts
Install .reg file	0	0
Create Key	10	10
Delete Key	0	0
Create/Set Value	25	25
Delete Value	0	0

Redistributable Packages	Attempts	Pass Counts
	0	0

Installation Qualification Result: PASS

Note: The installation qualification test passed because all installation actions were performed correctly (i.e., Attempts equals Pass Counts). The details of each file location, size, date and version, and registry setting are stored in file:

C:\Program Files (x86)\Thermo\IQ Reports\2022_04_07_115649 Thermo Exactive Series 2.12.0.3134 IQ Log.pdf

The following section shall be applied only to the optional qualification service. It is to be completed only by an individual that has received proper training in the installation of Thermo Fisher Scientific software, who affirms that this IQ Report is applicable only to the computer system that corresponds to that recorded in the respective software IQ worksheets for the Qualification number denoted on this page.

Trained Installer (Signature/Date): _____