




2695 PM Kit (WAT270944) Contents

Part	Part #	Qty.
Seal 1/8" Fluoroloy-G mini	WAT270789	2
Plunger	WAT270488	2
Solvent Filter	WAT25531	4
Check Valve: Cartridge Assy	289000144	2
Face Seal	WAT270163	4
Syringe	WAT073109	1
Filter Insert Assy	WAT088084	1
*Seal Wash Tube Seal	WAT270668	5
Seal Wash Plunger Seal Wash	WAT270160	2
Assy, Diffuser	WAT007272	1
* Screw, Compression	WAT025313	1
*Ferrule	410001125	1
*Spacer Upper Seal Wash	WAT031957	1
Battery, Lithium 484	WAT080443	1
*Needle Wash Frit	WAT072324	1
*Teflon Washer	WAT077321	1
*Tube, Cut Seal Wash Injector	WAT031956	1
*Needle	405001338	1
*Needle Seal	WAT033276	2
*Needle Wash Filter	415000230	1
*Retainer, Filter	405001171	1
Tubing, Teflon AWG 9 Thin Wall	WAT024036	1
*Glove, Powder Free, 2/pk	WAT241050	1

- The 2695 PM kit contains all the parts subject to normal wear to be replaced in one instrument.
- It comes with a 7 language protocol outlining instrument diagnostic tests, accuracy calibration, and tables for recording results.
- Items marked with an asterisk (*) are also found in the Seal Pack Rebuild Kit,
- P/N WAT271019.
- The kit contains one Seal Wash Tube Seal, P/N WAT270668.

