

Washer unit: PG 8593

Company: Waters Pacific Pte Ltd

Department: WPPL Chemical Prep Lab

Inventory no.: _____

Serial no.: 74400208

Table of Contents:

No.	Content
1	Documentation of Installation Qualification and Operation Qualification for the Miele Washer Unit PG 8593 and Peripheral Units
2	Certifications of equipment manufacturer
3	Documents/ certificates for washer unit
4	Certificate of technician
5	Test protocol of electrical safety test
6	Certificates of test equipment
7	Printout of settings of washer unit and programmes



Water Purification Systems





Water Purification Systems






sartorius

Equipment Qualification


Life Cycle Documentation



Qualification Protocol For Liebherr Products

Document Number: CQP-LIE-LKUEXV1610-01

REV: 1.0 Equipment Serial Number: 840303798



Qualification Protocol

Liebherr LKUEXV1610

CQP-LIE-LKUEXV1610-01

Version 1.0

Page 1 of 15

ITS SCIENCE & MEDICAL PTE. LTD.

Company Reg. No: 198700389G

219, Henderson Road, #11-02, Henderson Industrial Park, Singapore 159556.

Tel: +65-6273 0898 Fax: +65-6273 0810 E-mail: info@its-sciencemedical.com Website: <http://www.its-sciencemedical.com>


• Singapore • Malaysia • Thailand • Vietnam • Indonesia • Philippines • Myanmar

THE ITS GROUP

Qualification Protocol For Liebherr Products

Document Number: CQP-LIE-LGUEx1500-01

Rev: 1.0 Equipment Serial Number: 840442763



Qualification Protocol

Liebherr LGUEx1500

CQP-LIE-LGUEx1500-01

Version 1.0

Page 1 of 15

ITS SCIENCE & MEDICAL PTE. LTD.

Company Reg. No: 198700389G

219, Henderson Road, #11-02, Henderson Industrial Park, Singapore 159556.

Tel: +65-6273 0898 Fax: +65-6273 0810 E-mail: info@its-sciencemedical.com Website: <http://www.its-sciencemedical.com>

• Singapore • Malaysia • Thailand • Vietnam • Indonesia • Philippines • Myanmar

THE ITS GROUP

Performance Maintenance Protocol

515 Pump

Serial Number/ Firmware Version:	K06515 972A
Date:	11/OCT/2017
Performer:	KEVIN LAW

